

trans-9,10-Epoxy-octadecanoic acid, Purity: >98%, solution in Methanol , concentration: 10 mg/mL

Art. ID LA-12-1811-4
Unit 5 mg
Deliverydetails De Minimis UN1230 (6.1) II E2

Description

Synonyms: trans-9,10-Epoxystearic acid / (±)-trans-9,10-Epoxyoctadecanoic acid

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	trans-9,10-Epoxy-octadecanoic acid		10	mg/mL		