

Phosphatidylcholine, PC (egg) (in 1 ml CHCl₃), Purity: >98%, solution in Chloroform, concentration: Varies

Art. ID LA-37-0100-8S
Unit 50 mg
Deliverydetails De Minimis UN1888 III E1

Description

Synonyms: PC (egg) / L-?-phosphatidylcholine (egg)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Phosphatidylcholine, PC (egg) (in 1 ml CHCl ₃)	[97281-44-2]	Varies			