

**Phosphatidylcholine, PC (liver, bovine or porcine), Purity: >98%, solution in Chloroform, concentration: Varies**

Art. ID LA-37-0106-7  
Unit 25 mg  
Deliverydetails De Minimis UN1888 III E1

Description

Synonyms: PC (liver, bovine or porcine) / L-?-phosphatidylcholine (liver, bovine)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Phosphatidylcholine, PC (liver, bovine or porcine)	[475662-36-3]	Varies			