

Phosphatidylethanolamine, PE (corn germ), Purity: >98%, solution in Chloroform, concentration: 5 mg/mL

Art. ID LA-37-0122-4
Unit 5 mg
Deliverydetails De Minimis UN1888 III E1

Description

Synonyms: PE (corn germ) / L-?-phosphatidylethanolamine (corn germ)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Phosphatidylethanolamin e, PE (corn germ)		5	mg/mL		