

**Phosphatidylethanolamine, PE (liver, porcine), Purity: >98%, solution in Chloroform, concentration: Varies**

Art. ID LA-37-0126  
Unit on request  
Deliverydetails De Minimis UN1888 III E1

Description

Synonyms: PE (liver, porcine) / L-?-phosphatidylethanolamine (liver, porcine)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Phosphatidylethanolamin e, PE (liver, porcine)		Varies			