

Phosphatidylinositol, PI (soya), appr 50%, Purity: >98%, solution in Chloroform, concentration: Varies

Art. ID LA-37-0130-13
Unit 1 g
Deliverydetails De Minimis UN1888 III E1

Description

Synonyms: PI (plant) Ammonium salt / L-?-phosphatidylinositol (soybean)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Phosphatidylinositol, P I (soya), appr 50%		Varies			