

**Phosphatidic acid, PA (soybean) NH<sub>4</sub><sup>+</sup> salt, Purity: >98%, solution in Chloroform, concentration: Varies**

Art. ID LA-37-0141-7  
Unit 25 mg  
Deliverydetails De Minimis UN1888 III E1

Description

Synonyms: PA (soybean) NH<sub>4</sub><sup>+</sup> salt / L-?-phosphatidic acid (soybean) NH<sub>4</sub><sup>+</sup> salt

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Phosphatidic acid, PA (soybean) NH <sub>4</sub> <sup>+</sup> salt	[475995-54-1]	Varies			