

Phosphatidic acid, PA (soybean) NH₄⁺ salt, Purity: >98%, solution in Chloroform, concentration: Varies

Art. ID LA-37-0141-9
Unit 100 mg
Deliverydetails De Minimis UN1888 III E1

Description

Synonyms: PA (soybean) NH₄⁺ salt / L-?-phosphatidic acid (soybean) NH₄⁺ salt

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Phosphatidic acid, PA (soybean) NH ₄ ⁺ salt	[475995-54-1]	Varies			