

Phosphatidylserine, PS (brain, porcine), Purity: >98%, solution in Chloroform, concentration: Varies

Art. ID LA-37-0162-7
Unit 25 mg
Deliverydetails De Minimis UN1888 III E1

Description

Synonyms: PS (brain, porcine) / L-?-Phosphatidylserine (brain, bovine)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Phosphatidylserine, PS (brain, porcine)	[383907-32-2]	Varies			