

1,2-Diheptanoyl-sn-Glycerol-3-Phosphatidylcholine, PC(7:0/7:0), Purity: >98%

Art. ID LA-37-0700-0

Unit on request

Description

Synonyms: 7:0 PC / 1,2-diheptanoyl-sn-glycerol-3-phosphatidylcholine / [(2R)-2,3-di(heptanoyloxy)propyl] 2-trimethylazaniumylethyl phosphate

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	1,2-Diheptanoyl-sn-Glycerol-3-Phosphatidylcholine	[39036-04-9]				