

1,2-Dipalmitoyl-sn-Glycero-3-Phosphatidylcholine, PC(16:0/16:0), Purity: >99%

Art. ID LA-37-1600-13
Unit 1 g

Description

Synonyms: 16:0 PC / DPPC / 1,2-dihexadecanoyl-sn-glycero-3-phosphocholine /
[(2R)-2,3-di(hexadecanoyloxy)propyl] 2-trimethylazaniumylethyl phosphate

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	1,2-Dipalmitoyl-sn-Glycero-3-Phosphatidylcholine					