

**1,2-Diheptadecanoyl-sn-Glycero-3-Phosphatidylethanolamine, PE(17:0/17:0), Purity: >98%**

Art. ID LA-37-1720-0  
Unit on request

Description

Synonyms: 17:0 PE / 1,2-Diheptadecanoylphosphoethanolamine /  
1,2-diheptadecanoyl-sn-glycero-3-phosphoethanolamine

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	1,2-Diheptadecanoyl-sn- Glycero-3-Phosphatidyle thanolamine	[140219-78-9]				