

1,2-Distearoyl-sn-Glycerol-3-Phosphatidylcholine, PC(18:0/18:0), Purity: >98%

Art. ID LA-37-1800-7
Unit 25 mg

Description

Synonyms: 18:0 PC / DSPC / Distearoylphosphatidylcholine /
1,2-Bis(stearoyl)-sn-glycerol-3-phosphocholinylycerol-3-phosphocholine

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	1,2-Distearoyl-sn-Glyce ro-3-Phosphatidylcholin e					