

1,2-Distearoyl-sn-Glycerol-3-Phosphatidylethanolamine, PE(18:0/18:0), Purity: >99%

Art. ID LA-37-1820-7
Unit 25 mg

Description

Synonyms: 18:0 PE / DSPE / Distearoylphosphatidylethanolamine /
1,2-dioctadecyl-rac-glycerol-3-phosphoethanolamine

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	1,2-Distearoyl-sn-Glyce ro-3-Phosphatidylethano lamine					