

1,2-Diheptadecanoin-sn-Glycero-3-Phosphatidylserine Na salt, PS(17:0/17:0), Purity: >98%

Art. ID LA-38-1780-0

Unit on request

Description

Synonyms: 17:0 PS / Sodium 1,2-diheptadecanoyl-sn-glycero-3-PS / 1,2-diheptadecanoyl-sn-glycero-3-phospho-L-serine (sodium salt)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	1,2-Diheptadecanoin-sn-Glycero-3-Phosphatidylserine Na salt	[799268-51-2]				