

**1,2-Distearoyl-sn-Glycerol-3-Phosphatidic acid Na salt, PA(18:0/18:0), Purity: >99%**

Art. ID LA-38-1860-11  
Unit 250 mg

Description

Synonyms: 18:0 PA / DSPA / 1,2-distearoyl-sn-glycerol-3-phosphatidic acid (sodium salt) /  
1,2-dioctadecanoyl-sn-glycerol-3-phosphate (sodium salt)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	1,2-Distearoyl-sn-Glyce ro-3-Phosphatidic acid Na salt					