

**Single Crystal Diamond Probes specially designed for Nanoindentation, Nanoscratching, Lithography, resonant frequency 400-650kHz, force constant 100-200 N/m, curvature radius <15nm**

Art. ID	TIPS-DRP_IN
Unit	Set of 5 probes
Deliverydetails	No Dangerous Good /not restricted

Description

Unique probes with Single Crystal Diamond tip specially designed for Nanoindentation, Nanoscratching, Lithography etc. Resonant frequency 400-650kHz, force constant 100-200 N/m, curvature radius <15nm.