

**TOP VISUAL High Resolution Contact Silicon AFM Cantilevers VIT\_P\_C-A series are specially designed for tip or laser spot precise positioning over the point of interest**

Art. ID	TIPS-VIT_P_C-A/50
Unit	Set of 50 cantilevers
Deliverydetails	No Dangerous Good /not restricted

Description

TOP VISUAL High Resolution Contact Silicon AFM Cantilevers VIT\_P\_C-A series are specially designed for tip or laser spot precise positioning over the point of interest. Typical Resonant Frequency 16 kHz (guaranteed range 8-25 kHz), Typical Force Constant 0.3 N/m (guaranteed range 0.06-1N/m). Cantilevers have Al reflective coating.