

**High Resolution Semicontact / Noncontact Silicon AFM Cantilevers NSG01 series,
resonant frequency 87-230kHz, force constant 1,45-15,1N/m - Separated chips**

Art. ID	TIPS-NSG01/Pt/200
Unit	200 separated chips
Deliverydetails	No Dangerous Good /not restricted

Description

High Resolution Silicon AFM Cantilevers NSG01 series are designed for Semicontact (Intermittent), Noncontact applications. Typical Resonant Frequency 150 kHz (guaranteed range 87-230kHz), Typical Force Constant 5.1 N/m (guaranteed range 1.45-15.1N/m). Cantilever has Au reflective side coating to increase laser signal.