

**High Resolution Semicontact / Noncontact Silicon AFM Cantilevers NSG30 series, resonant frequency 240-440 kHz, force constant 22-100 N/m - Separated chips**

Art. ID	TIPS-NSG30/Pt/200
Unit	200 separated chips
Deliverydetails	No Dangerous Good /not restricted

Description

High Resolution Silicon AFM Cantilevers NSG30 series are designed for Semicontact ( Intermittent ), Noncontact applications. Typical Resonant Frequency 320 kHz (guaranteed range 240-440 kHz), Typical Force Constant 40 N/m (guaranteed range 22-100 N/m). Cantilever has Au reflective side coating to increase laser signal.