

**High Resolution Force Modulation Silicon AFM Cantilevers FMG01 series, resonant frequency 40-96kHz, force constant 1-3N/m - Separated chips**

Art. ID	TIPS-FMG01/Au/100
Unit	100 separated chips
Deliverydetails	No Dangerous Good /not restricted

Description

High Resolution Silicon AFM Cantilevers FMG01 series are designed for Force Modulation applications. Typical Resonant Frequency 60 kHz (guaranteed range 40-96kHz), Typical Force Constant 3 N/m (guaranteed range 1-5N/m). Cantilever has Au reflective side coating to increase laser signal.