

**High Resolution Contact Silicon AFM Cantilevers CSG01 series, resonant frequency 4-17kHz, force constant 0,003-0,13N/m - Separated chips**

Art. ID	TIPS-CSG01/Pt/50
Unit	50 separated chips
Deliverydetails	No Dangerous Good /not restricted

Description

High Resolution Silicon AFM Cantilevers CSG01 series are designed for Contact applications. Typical Resonant Frequency 9,8 kHz (guaranteed range 4-17kHz), Typical Force Constant 0,03 N/m (guaranteed range 0,003-0,13N/m). Cantilever has Au reflective side coating to increase laser signal.