

**High Resolution Semicontact / Noncontact Silicon AFM Cantilevers NSG03 series, resonant frequency 47-150kHz, force constant 0.35-6.1N/m - Separated chips**

Art. ID	TIPS-NSG03/Pt/50
Unit	50 separated chips
Deliverydetails	No Dangerous Good /not restricted

Description

High Resolution Silicon AFM Cantilevers NSG03 series are designed for Semicontact ( Intermittent ), Noncontact applications. Typical Resonant Frequency 90 kHz (guaranteed range 47-150kHz), Typical Force Constant 1.74 N/m (guaranteed range 0.35-6.1N/m). Cantilever has Au reflective side coating to increase laser signal.