

High Resolution Contact Silicon AFM Cantilevers CSG10 series, resonant frequency 8-39kHz, force constant 0,01-0,5N/m - Separated chips

Art. ID	TIPS-CSG10/Pt/15
Unit	15 separated chips
Deliverydetails	No Dangerous Good /not restricted

Description

High Resolution Silicon AFM Cantilevers CSG10 series are designed for Contact applications. Typical Resonant Frequency 22 kHz (guaranteed range 8-39kHz), Typical Force Constant 0,11 N/m (guaranteed range 0,01-0,5N/m). Cantilever has Au reflective side coating to increase laser signal