

Diamond Coated Semicontact / Noncontact Conductive Probes DCP01 series, resonant frequency 115-190 kHz, force constant 2,5-10 N/m

Art. ID	TIPS-DCP01/50
Unit	50 separated chips
Deliverydetails	No Dangerous Good /not restricted

Description

Diamond Like Coated Conductive Probes DCP01 series are The best choice for any kind of long-term electrical characterizations. Stable and nondestructive, wear resistant probes with conductive diamond-like coating allow you to make as many images as you want. Resonant frequency 115-190 kHz, force constant 2,5-10 N/m. Cantilever has Au reflective side coating to increase laser signal