

Diamond Coated Semicontact / Noncontact Conductive Probes DCP10 series, resonant frequency 190-325 kHz, force constant 5,5-22.5 N/m

Art. ID	TIPS-DCP10/15
Unit	15 separated chips
Deliverydetails	No Dangerous Good /not restricted

Description

Diamond Coated Conductive Probes DCP10 series are the best choice for any kind of long-term electrical characterizations. Stable and nondestructive, wear resistant probes with conductive diamond-like coating allow you to make as many images as you want. Resonant frequency 190-325 kHz, force constant 5,5-22.5 N/m. Cantilever has Au reflective side coating to increase laser signal