

**High Resolution High Accuracy Semicontact / Noncontact / Force Modulation AFM Cantilevers HA FM series with Pt tip conductive coating, each chip has 2 cantilevers, resonant frequency 114 kHz / 77 kHz, force constant 6 N/m / 3.5 N/m**

|                 |                                   |
|-----------------|-----------------------------------|
| Art. ID         | TIPS-HA_FM/Pt/300                 |
| Unit            | 300 separated chips               |
| Deliverydetails | No Dangerous Good /not restricted |