

## Cantilever set

Art. ID	TIPS-UNISET/25
Unit	set
Deliverydetails	No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Cantilever set consists						
of						
5 probes for noncontact						
/semicontact modes - N						
SG30						
5 probes for contact mo						
des - CSG10						
5 probes for PFM - FMG0						
1/Pt						
5 probes for SRM - CSG1						
0/Au						
5 probes for EFM, SCM,						
SKM - DCP10						