

Cantilever set

Art. ID	TIPS-UNISET/50
Unit	set
Deliverydetails	No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Cantilever set consists of						
10 probes for noncontact/semicontact modes - N SG10						
10 probes for force modulation mode - FMG01						
10 probes for contact mode - CSG10						
5 probes for EFM, SCM, SKM modes - NSG01/Pt						
5 probes for SRM - CSG10/Pt						
5 probes with stable and nondestructive diamond like conductive coating for LAO Lithography and other conductive modes - DCP20						