

## **AFM probe sets**

Art. ID	TIPS-SEMISET/47
Unit	set
Deliverydetails	No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
AFM probe sets:						
10 chips - for contact						
AFM						
20 chips - for noncontact modes						
5 chips - for spreading resistance						
5 chips - for SKM/SCM						
5 chips - for MFM						
2 chips NSC05 -Whisker type- for deep trenches measurements						