

## **ARCONIC - Hydrogen in Aluminium Reference Material, rod, 5inch x 24inch**

Art. ID ALCOA-Alu-H2  
Unit 10 rods  
Deliverydetails No Dangerous Good /not restricted

### Description

Certified: January 15, 2014 /// ALCOA-Alu-H2 is provided in rod form and is intended for use in inert gas fusion analysis of hydrogen in aluminium as described in ASTM method E2792.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Hydrogen (H)	[1333-74-0]	0,068	mg/kg	ASTM E2792	