

Custom Solution - 12 components in 5% HNO₃

Art. ID LM24-CUS-70166
Unit 2 x 250 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Conservation 24 months (12 months after opening)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Vanadium (V)	[7440-62-2]	1000	mg/L		
	Arsenic (As)	[7440-38-2]	1000	mg/L		
	Cadmium (Cd)	[7440-43-9]	1000	mg/L		
	Cobalt (Co)	[7440-48-4]	1000	mg/L		
	Copper (Cu)	[7440-50-8]	1000	mg/L		
	Lead (Pb)	[7439-92-1]	1000	mg/L		
	Zinc (Zn)	[7440-66-6]	1000	mg/L		
	Selenium (Se)	[7782-49-2]	1000	mg/L		
	Nickel (Ni)	[7440-02-0]	1000	mg/L		
	Thallium (Tl)	[7440-28-0]	1000	mg/L		
	Strontium (Sr)	[7440-24-6]	1000	mg/L		
	Silver (Ag)	[7440-22-4]	1000	mg/L		