

### **Custom Solution - 5 components in 2% HNO<sub>3</sub>**

Art. ID LM24-CUS-70168  
Unit 2 x 100 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

#### **Description**

Conservation 24 months (12 months after opening)

| Text/Information | Analyte/Parameter | CAS number  | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
|                  | Molybdenum (Mo)   | [7439-98-7] | 1000                | mg/L |        |        |
|                  | Zirconium (Zr)    | [7440-67-7] | 1000                | mg/L |        |        |
|                  | Tungsten (W)      | [7440-33-7] | 1000                | mg/L |        |        |
|                  | Hafnium (Hf)      | [7440-58-6] | 1000                | mg/L |        |        |
|                  | Niobium (Nb)      | [7440-03-1] | 1000                | mg/L |        |        |