

Custom Solution - 5 components in 2% HNO₃

Art. ID LM24-CUS-70169
Unit 2 x 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Description

Conservation 12 months (6 months after opening)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Molybdenum (Mo)	[7439-98-7]	10	mg/L		
	Zirconium (Zr)	[7440-67-7]	10	mg/L		
	Tungsten (W)	[7440-33-7]	10	mg/L		
	Hafnium (Hf)	[7440-58-6]	10	mg/L		
	Niobium (Nb)	[7440-03-1]	10	mg/L		