

Iron Base, Cast Iron, RM, Mushroom 2 pieces 43 mm x 5 mm

Art. ID CTIF-FL10
Unit 2 pieces
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	~0,03	%		
	Arsenic (As)	[7440-38-2]	~0,022	%		
	Soluble B		~0,012	%		
	Bismuth (Bi)	[7440-69-9]	~0,004	%		
	Carbon (C)	[7440-44-0]	3,1	%		
	Chromium (Cr)	[7440-47-3]	~0,07	%		
	Copper (Cu)	[7440-50-8]	0,104	%		
	Manganese (Mn)	[7439-96-5]	0,85	%		
	Molybdenum (Mo)	[7439-98-7]	0,0335	%		
	Niobium (Nb)	[7440-03-1]	~0,018	%		
	Nickel (Ni)	[7440-02-0]	0,1	%		
	Phosphorus (P)	[7723-14-0]	0,323	%		
	Lead (Pb)	[7439-92-1]	~0,002	%		
	Sulfur (S)	[7704-34-9]	0,066	%		
	Antimony (Sb)	[7440-36-0]	~0,032	%		
	Silicon (Si)	[7440-21-3]	1,3	%		
	Tin (Sn)	[7440-31-5]	0,028	%		
	Tellurium (Te)	[13494-80-9]	~0,001	%		
	Titanium (Ti)	[7440-32-6]	0,045	%		
	Vanadium (V)	[7440-62-2]	0,048	%		
	Tungsten (W)	[7440-33-7]	~0,02	%		
	Zinc (Zn)	[7440-66-6]	~0,029	%		