

## **CONOSTAN - Multi-element standard in 75 blank oil: S21 + K:500 ppm**

Art. ID                      CON-150-021-059  
Unit                         200 g  
Deliverydetails            No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number  | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
|                  | Aluminium (Al)    | [7429-90-5] | 0,05                | %    |        |        |
|                  | Aluminium (Al)    | [7429-90-5] | 500                 | µg/g |        |        |
|                  | Barium (Ba)       | [7440-39-3] | 0,05                | %    |        |        |
|                  | Barium (Ba)       | [7440-39-3] | 500                 | µg/g |        |        |
|                  | Boron (B)         | [7440-42-8] | 0,05                | %    |        |        |
|                  | Boron (B)         | [7440-42-8] | 500                 | µg/g |        |        |
|                  | Cadmium (Cd)      | [7440-43-9] | 0,05                | %    |        |        |
|                  | Cadmium (Cd)      | [7440-43-9] | 500                 | µg/g |        |        |
|                  | Calcium (Ca)      | [7440-70-2] | 0,05                | %    |        |        |
|                  | Calcium (Ca)      | [7440-70-2] | 500                 | µg/g |        |        |
|                  | Chromium (Cr)     | [7440-47-3] | 0,05                | %    |        |        |
|                  | Chromium (Cr)     | [7440-47-3] | 500                 | µg/g |        |        |
|                  | Copper (Cu)       | [7440-50-8] | 0,05                | %    |        |        |
|                  | Copper (Cu)       | [7440-50-8] | 500                 | µg/g |        |        |
|                  | Iron (Fe)         | [7439-89-6] | 0,05                | %    |        |        |
|                  | Iron (Fe)         | [7439-89-6] | 500                 | µg/g |        |        |
|                  | Lead (Pb)         | [7439-92-1] | 0,05                | %    |        |        |
|                  | Lead (Pb)         | [7439-92-1] | 500                 | µg/g |        |        |
|                  | Magnesium (Mg)    | [7439-95-4] | 0,05                | %    |        |        |
|                  | Magnesium (Mg)    | [7439-95-4] | 500                 | µg/g |        |        |
|                  | Manganese (Mn)    | [7439-96-5] | 0,05                | %    |        |        |
|                  | Manganese (Mn)    | [7439-96-5] | 500                 | µg/g |        |        |
|                  | Molybdenum (Mo)   | [7439-98-7] | 0,05                | %    |        |        |
|                  | Molybdenum (Mo)   | [7439-98-7] | 500                 | µg/g |        |        |
|                  | Nickel (Ni)       | [7440-02-0] | 0,05                | %    |        |        |
|                  | Nickel (Ni)       | [7440-02-0] | 500                 | µg/g |        |        |
|                  | Phosphorus (P)    | [7723-14-0] | 0,05                | %    |        |        |
|                  | Phosphorus (P)    | [7723-14-0] | 500                 | µg/g |        |        |
|                  | Potassium (K)     | [7440-09-7] | 0,05                | %    |        |        |
|                  | Potassium (K)     | [7440-09-7] | 500                 | µg/g |        |        |
|                  | Silicon (Si)      | [7440-21-3] | 0,05                | %    |        |        |

|               |             |      |      |
|---------------|-------------|------|------|
| Silicon (Si)  | [7440-21-3] | 500  | µg/g |
| Silver (Ag)   | [7440-22-4] | 0,05 | %    |
| Silver (Ag)   | [7440-22-4] | 500  | µg/g |
| Sodium (Na)   | [7440-23-5] | 0,05 | %    |
| Sodium (Na)   | [7440-23-5] | 500  | µg/g |
| Tin (Sn)      | [7440-31-5] | 0,05 | %    |
| Tin (Sn)      | [7440-31-5] | 500  | µg/g |
| Titanium (Ti) | [7440-32-6] | 0,05 | %    |
| Titanium (Ti) | [7440-32-6] | 500  | µg/g |
| Vanadium (V)  | [7440-62-2] | 0,05 | %    |
| Vanadium (V)  | [7440-62-2] | 500  | µg/g |
| Zinc (Zn)     | [7440-66-6] | 0,05 | %    |
| Zinc (Zn)     | [7440-66-6] | 500  | µg/g |