

Art. ID GBW-01403a
Unit

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Aluminium (Al) | [7429-90-5] | 0,06 | % | | |
| | Carbon (C) | [7440-44-0] | 0,012 | % | | |
| | Chromium (Cr) | [7440-47-3] | 0,016 | % | | |
| | Copper (Cu) | [7440-50-8] | 0,021 | % | | |
| | Manganese (Mn) | [7439-96-5] | 0,034 | % | | |
| | Nitrogen (N) | [7727-37-9] | 0,015 | % | | |
| | Nickel (Ni) | [7440-02-0] | 0,03 | % | | |
| | Phosphorus (P) | [7723-14-0] | 0,0019 | % | | |
| | Sulfur (S) | [7704-34-9] | 0,02 | % | | |
| | Silicon (Si) | [7440-21-3] | 0,017 | % | | |