

Cannabichromene-D3 (CBC-D3), 100 µg/mL in methanol, certified reference material, ampule of 1 mL

Art. ID SAF-C-211-1ML
Unit 1 x 1 ml
Deliverydetails De Minimis UN1230 (6.1) II E2

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-----------------------------|------------|---------------------|-------|--------|--------|
| | Cannabichromene-D3 (CBC-D3) | | 100 | µg/mL | | |