

Molindone-D8 HCl, 100 µg/mL in methanol (as free base), certified reference material, ampule of 1 m

Art. ID SAF-M-235-1ML
Unit 1 x 1 ml
Deliverydetails De Minimis UN1230 (6.1) II E2

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|----------------|---------------------|-------|--------|--------|
| | Molindone-D8 HCl | [1189669-25-7] | 100 | µg/mL | | |