

### **Custom Solution - 7 components in HNO3 5%, HF 0.8%**

Art. ID LM24-CUS-71295  
Unit 125 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN2922 (6.1) II E2

#### Description

Conservation: 12 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silicon (Si)	[7440-21-3]	1000	mg/L		
	Rhenium (Re)	[7440-15-5]	100	mg/L		
	Phosphorus (P)	[7723-14-0]	500	mg/L		
	Antimony (Sb)	[7440-36-0]	500	mg/L		
	Tin (Sn)	[7440-31-5]	500	mg/L		
	Tellurium (Te)	[13494-80-9]	500	mg/L		
	Molybdenum (Mo)	[7439-98-7]	500	mg/L		