

### **Custom Solution - 3 components in HNO3 5%**

Art. ID LM24-CUS-71300  
Unit 1 x 100 mL bottle  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

#### Description

Conservation 24 months. (12 months after opening)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Tin (Sn)	[7440-31-5]	5455	µg/mL		
	Molybdenum (Mo)	[7439-98-7]	2727	µg/mL		
	Antimony (Sb)	[7440-36-0]	1091	µg/mL		