

Custom Solution - 3 components in HNO3 5%

Art. ID LM24-CUS-71301
Unit 1 x 250 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Conservation 24 months. (12 months after opening)

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|-------|--------|--------|
| | Tin (Sn) | [7440-31-5] | 5455 | µg/mL | | |
| | Molybdenum (Mo) | [7439-98-7] | 2727 | µg/mL | | |
| | Antimony (Sb) | [7440-36-0] | 1091 | µg/mL | | |