

## **Custom Solution - 7 components in HNO3 10%**

Art. ID LM24-CUS-71303  
Unit 1 x 100 mL bottle  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN2031 II E2

### Description

Conservation 12 months. (6 months after opening)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Potassium (K)	[7440-09-7]	500	mg/L		
	Nickel (Ni)	[7440-02-0]	500	mg/L		
	Copper (Cu)	[7440-50-8]	100	mg/L		
	Manganese (Mn)	[7439-96-5]	100	mg/L		
	Barium (Ba)	[7440-39-3]	20	mg/L		
	Calcium (Ca)	[7440-70-2]	20	mg/L		
	Zinc (Zn)	[7440-66-6]	20	mg/L		