

## **Custom Solution - 7 components in HNO3 2%**

Art. ID LM24-CUS-71310  
Unit 100 mL bottle  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

### Description

Conservation 24 months. (12 months after opening)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Antimony (Sb)	[7440-36-0]	100	mg/L		
	Molybdenum (Mo)	[7439-98-7]	100	mg/L		
	Silicon (Si)	[7440-21-3]	1000	mg/L		
	Titanium (Ti)	[7440-32-6]	100	mg/L		
	Boron (B)	[7440-42-8]	100	mg/L		
	Phosphorus (P)	[7723-14-0]	1000	mg/L		
	Tin (Sn)	[7440-31-5]	100	mg/L		