

### **Custom Solution - 9 components 1000 mg/L each of in HNO3 5%**

Art. ID LM24-CUS-71328  
Unit 1 x 100 mL bottle  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

#### **Description**

Conservation 24 months (12 months after opening)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Antimony (Sb)	[7440-36-0]	1000	mg/L		
	Arsenic (As)	[7440-38-2]	1000	mg/L		
	Beryllium (Be)	[7440-41-7]	1000	mg/L		
	Molybdenum (Mo)	[7439-98-7]	1000	mg/L		
	Selenium (Se)	[7782-49-2]	1000	mg/L		
	Tin (Sn)	[7440-31-5]	1000	mg/L		
	Tellurium (Te)	[13494-80-9]	1000	mg/L		
	Vanadium (V)	[7440-62-2]	1000	mg/L		
	Titanium (Ti)	[7440-32-6]	1000	mg/L		