

Canola BASF Event T45 Leaf DNA 0208-A7

Art. ID	AOCS-0208-A7
Unit	approx. 10 µg
Deliverydetails	No Dangerous Good /not restricted

Description

AOCS-0208-A7 was prepared from canola genomic DNA extracted from clean T45 canola leaves provided by BASF Agricultural Solutions Seed US LLC. AOCS-0208-A7 is available in sterile, 0.5 ml skirted, screw-cap, self-sealing tubes containing approximately 10 µg extracted DNA.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
T45 canola leaf tissue genomic DNA	T45 DNA		>999,99	ng/µg	event-spec ific real- time PCR	