

## **Custom Solution - 21 components in HNO3 2%**

Art. ID LM24-CUS-71423  
Unit 1 x 100 mL bottle  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

### Description

Conservation 12 months (originally closed bottle). (6 months after opening)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	200	mg/L		
	Arsenic (As)	[7440-38-2]	100	mg/L		
	Boron (B)	[7440-42-8]	20	mg/L		
	Barium (Ba)	[7440-39-3]	200	mg/L		
	Calcium (Ca)	[7440-70-2]	20000	mg/L		
	Cadmium (Cd)	[7440-43-9]	20	mg/L		
	Cobalt (Co)	[7440-48-4]	100	mg/L		
	Chromium (Cr)	[7440-47-3]	100	mg/L		
	Copper (Cu)	[7440-50-8]	100	mg/L		
	Iron (Fe)	[7439-89-6]	2000	mg/L		
	Potassium (K)	[7440-09-7]	10000	mg/L		
	Magnesium (Mg)	[7439-95-4]	10000	mg/L		
	Manganese (Mn)	[7439-96-5]	100	mg/L		
	Molybdenum (Mo)	[7439-98-7]	100	mg/L		
	Sodium (Na)	[7440-23-5]	10000	mg/L		
	Nickel (Ni)	[7440-02-0]	100	mg/L		
	Phosphorus (P)	[7723-14-0]	2000	mg/L		
	Thallium (Tl)	[7440-28-0]	200	mg/L		
	Vanadium (V)	[7440-62-2]	100	mg/L		
	Zinc (Zn)	[7440-66-6]	200	mg/L		
	Lead (Pb)	[7439-92-1]	100	mg/L		