

**Custom Solution - 1 component, TI 24 mg/L in HNO3 2%**

Art. ID LM24-CUS-71426  
Unit 1 x 125 mL bottle  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

**Description**

Conservation 12 months (originally closed bottle). (6 months after opening)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Thallium (TI)	[7440-28-0]	24	mg/L		