

PAH Mix, 18 compounds., 100 µg/mL each of in in Dichloromethane / Methanol (1/1)

Art. ID LM24-OP-4011667
Unit 1 x 1 mL ampoule
Deliverydetails De Minimis

Description

Conservation: 24 months (originally closed bottle).

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|------------------------|------------|---------------------|-------|--------|--------|
| | Benzo(b)fluoranthene | [205-99-2] | 100 | µg/mL | | |
| | Chrysene | [218-01-9] | 100 | µg/mL | | |
| | Dibenz(a,h)anthracene | [53-70-3] | 100 | µg/mL | | |
| | Fluoranthene | [206-44-0] | 100 | µg/mL | | |
| | Fluorene | [86-73-7] | 100 | µg/mL | | |
| | Indeno(1,2,3-cd)pyrene | [193-39-5] | 100 | µg/mL | | |
| | Naphthalene | [91-20-3] | 100 | µg/mL | | |
| | Pyrene | [129-00-0] | 100 | µg/mL | | |
| | Phenanthrene | [85-01-8] | 100 | µg/mL | | |
| | Benzo(e)pyrene | [192-97-2] | 100 | µg/mL | | |
| | Perylene | [198-55-0] | 100 | µg/mL | | |
| | Acenaphthene | [83-32-9] | 100 | µg/mL | | |
| | Acenaphthylene | [208-96-8] | 100 | µg/mL | | |
| | Anthracene | [120-12-7] | 100 | µg/mL | | |
| | Benzo(a)pyrene | [50-32-8] | 100 | µg/mL | | |
| | Benzo(g,h,i)perylene | [191-24-2] | 100 | µg/mL | | |
| | Benzo(k)fluoranthene | [207-08-9] | 100 | µg/mL | | |
| | Benzo(a)anthracene | [56-55-3] | 100 | µg/mL | | |