

USP <541> Potentiometric Titration Proficiency Test

Art. ID USP-PHPT-004

Unit each

Description

Quantitative Proficiency Test (PT) for laboratories performing Potentiometric Titration Assay by USP <541> or equivalent. Participants report quantitative results for their sample(s). Samples provided two vials containing 500mg of material. Participants are asked to determine the concentration (Assay) in percentage (%) of target analyte in the sample per the USP monograph or equivalent. /// Target Analyte: Sodium Chloride, powdered or granulated solid

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Target Analyte, powdered or granulated solid	Sodium Chloride	[7647-14-5]				